Uro-Tainer® M

NaCl 0,9% B. Braun

BG	Указания за употреба	2
CZ	Informace pro použití	4
DE	Gebrauchsanweisung	6
DK	Brugsanvisning	8
ES	Instrucciones de uso	10
FI	Käyttöohje	12
FR	Mode d'emploi	14
GB	Information for use	16
HU	Használati utasítás	18
IT	Istruzioni per l'uso	20
NL	Gebruiksaanwijzing	22
NO	Brukerveiledning	24
PT	Instruções de uso	26
SE	Bruksanvisning	28
SK	Informácie pre použitie	30

schwarz

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Lätus - - -

Irrigation Solution__190_GIF
190/NU19091/1212
GIF (A5-Booklet, 15-spr.)
Standort Crissier

B | BRAUN

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Description

Uro-Tainer® NaCl 0.9% «M»

Solution for intravesical instillation of medication such as antibiotics, disinfectants, cytotoxics, antimycotics, corticosteroids, etc, via an indwelling or single-use urinary catheter.

Uro-Tainer® is a ready-to-use, disposable system. It consists of a PVC sachet, in a protective plastic overwrapping, containing 50 ml or 100 ml of sterile solution, with a separate port for adding medication, a flexible tube fitted with a clamp, and a sterile universal connector protected by a rigid plastic cap with a tamperproof seal.

Dosage

As required.

Overdosage

No clinical data exists concerning overdosage.

Important information

- This solution is intended for urinary catheter rinsing only : not for intravenous infusion.
- Instil the solution by gravity feed, avoid force.
- Precautions should be taken when performing bladder washout on patients with spinal injuries because of the risk of autonomic dysreflexia.
- Apply using aseptic technique. A port disinfection before use is recommended.
- When adding medication, refer to the manufacturer's instructions and verify its compatibility with NaCl 0.9%.
- When adding cytotoxics, it is essential to refer to the manufacturer's instructions special precautions: contraindications, preparation, administration, etc.
- Only to be used if the contents and the sealing cap are intact, and the solution is clear.
- This is a single-use system. After opening, it must be used and disposed of promptly.

Uro-Tainer® M

Maintenance solution for urethral and suprapubic catheters

NaCl 0,9% B. Braun

Sterile, free of endotoxins

- Dispose of any unused solution and the sachet. If medication has been added, refer to the manufacturer's disposal instructions.
- Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.
- Animal experiments have shown DEHP to be potentially toxic to reproduction. Given the current state of scientific knowledge, a risk (in case of long-term use) especially for premature male babies cannot be completely excluded. As a precaution, the application of medical products containing DEHP should be restricted to short-term use for pregnant women, breastfeeding mothers, infants, and children.
- Store out of the reach of children.
- Refer to the illustrated instructions for use on the next page of this leaflet.
- Discard properly after use.

Secondary effects

This solution is not known to have any secondary effects. If the liquid exits from the urethra during the irrigation, the catheter may no longer be in the bladder and the patient may require a new catheter.

If the liquid does not flow, the catheter may be kinked in the bladder or it may be blocked and may require changing.

Expiry date

Do not use after the date marked EXP on the back of the sachet.

Storage


Store in the original packaging, away from light and at a temperature not exceeding 25 °C.

Latest Information update


06/2010

The Uro-Tainer® M sachet is supplied sterile and should not be removed from its overwrapping until required for use. If desired, bring the solution to body temperature by immersing the wrapped sachet in warm water.


1. Tear open plastic bag and take out Uro-Tainer® M


2. Ensure Uro-Tainer® M tube is closed by plastic clamp.


3a. Using aseptic technique, assemble transfer needle and vial onto the injection port. Squeeze saline into the vial from the Uro-Tainer® M. Release pressure from the Uro-Tainer® M to draw the contents of the vial into the sachet. The squeezing and release of air from the Uro-Tainer® M may be repeated several times until all the solution is safely transferred; or


3b. Syringe some saline out of the Uro-Tainer® M through the injection port and into the vial to reconstitute the drug. Draw the reconstituted drug back into the syringe and inject into the Uro-Tainer® M again using the injection port. If a liquid drug is used, syringe it directly into Uro-Tainer® M through the injection port.

4. Remove tamper proof seal. Twist and withdraw cap without touching sterile connector.


5. Insert connector into the catheter funnel.

6. Holding Uro-Tainer® M above the level of the bladder, release the clamp. Apply the solution by gravity alone.


7. Close tube with clamp of the Uro-Tainer® M.

8. After the required period, release clamp and hold bag to have below the level of the bladder, drain the solution. Uro-Tainer® M can contain up to 150 ml and therefore allows to collect not only the irrigation solution but also eventual urine from the bladder.

9. Clamp tube and remove the connector from catheter by turning gently. Connect catheter to a sterile drainage bag.

Do not re-use, for single use only. Discard after treatment.

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B. Braun Medical AG